

## ClearLLab 10C Panels for leukemia and lymphoma(1st\_line)

Tube	FITC	PE	ECD	PC5.5	PC7	APC	APC-A700	APC-A750	PB	KRO
B Cell Tube	Kappa	Lambda	CD10	CD5	CD200	CD34	CD38	CD20	CD19	CD45
T Cell Tube	TCR $\gamma\delta$	CD4	CD2	CD56	CD5	CD34	CD7	CD8	CD3	CD45
M1 Cell Tube	CD16	CD7	CD10	CD13	CD64	CD34	CD14	HLA-DR	CD11b	CD45
M2 Cell Tube	CD15	CD123	CD117	CD13	CD33	CD34	CD38	HLA-DR	CD19	CD45

### 10色抗體管套組 ClearLLab 10C panel(第一線)

ClearLLab 10C B: kappa/lambda/CD10/CD5/CD200/CD34/CD38/CD20/CD19/CD45

ClearLLab 10C T: TCR $\gamma\delta$ /CD4/CD2/CD56/CD5/CD34/CD7/CD8/CD3/CD45

ClearLLab 10C M1: CD16/CD7/CD10/CD13/CD64/CD34/CD14/HLA-DR/CD11b/CD45

ClearLLab 10C M2: CD15/CD123/CD117/CD13/CD33/CD34/CD38/HLA-DR/CD11b/CD45

## LUCID DURA Innovations Panels for leukemia and lymphoma(2nd\_line)

	FITC	PE	ECD	PC5.5	PC7	APC	APC-A750	PB	KRO
ALOT	cyMPO	cyCD79a	CD34		CD19	CD7	SmCD3	CyCD3	CD45
B-CLPD-1	CD23	CD10		CD79b	CD19	CD200	CD43	CD20	CD45
MM-1	CD38	CD56		$\beta$ 2micro	CD19	cykappa	cylambda	CD138	CD45
MM-2	CD38	CD28		CD27	CD19	CD117	CD81	CD138	CD45

### 8色抗體管套組 DURA Innovations panel(第二線·依照疾病類別選擇性加染)

ALOT: cyMPO/cyCD79a/CD34/CD19/CD7/SmCD3/CyCD3/CD45

B-CLPD-1: CD23/CD10/CD79b/CD19/CD200/CD43/CD20/CD45

MM1: CD38/CD56/ $\beta$ 2 microglobulin /CD19/CyKappa/CyLambda/CD138/CD45

MM2: CD38/CD28/CD27/CD19/CD117/CD81/CD138/CD45